Guía Referencia Protocolo RM

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GE 1.5T MR Exam Protocol

Name: 3-plane Localizer

Pulse Sequence: 2D GRE
 Scan plane: axial, coronal, sagittal
 Slices: 5 coronal, 5 sagittal, 7 axial (17 total)
 TR: 60msec
 TE: 10msec
 Flip angle: 50°
 TI: N/A
 Echo train length: N/A
 NEX: 1
 Slice thickness: 5mm
 Skip: 1mm
 Matrix: 256 x 128
 FOV: 24cm
 Bandwidth: +/- 32 kHz
 Readout direction: Sagittal S/I, Coronal S/I, Axial A/P
 Options: none

Name: Coronal 2D PD FS

Pulse Sequence: 2D FSE
 Scan plane: coronal
 Slices: 28
 TR: 4500msec
 TE: Minimum full TE (~15msec)
 Flip angle: N/A
 TI: N/A
 Echo train length: 8
 NEX: 2
 Slice thickness: 3mm
 Skip: 0.5mm
 Matrix: 512 x 256
 FOV: 14cm
 Bandwidth: +/- 20-22 kHz
 Readout direction: S/I
 Options: Fat Saturation, no phase wrap, fast, TRF, EDR, ZIP512

Name: Sagittal 2D PD FS

Pulse Sequence: 2D FSE
 Scan plane: sagittal
 Slices: 28
 TR: 4500msec
 TE: Minimum full TE (~15msec)
 Flip angle: N/A
 TI: N/A
 Echo train length: 8
 NEX: 2
 Slice thickness: 3mm
 Skip: 0.5mm
 Matrix: 512 x 256
 FOV: 14cm
 Bandwidth: +/- 20-22 kHz
 Readout direction: S/I
 Options: Fat Saturation, no phase wrap, fast, TRF, EDR, ZIP512

Name: Sagittal 3D SPGR FS

Pulse Sequence: 3D SPGR
 Scan plane: sagittal
 Slices: 32
 TR: 60msec
 TE: 5msec
 Flip angle: 40°
 TI: N/A
 Echo train length: N/A
 NEX: 1
 Slice thickness: 3mm
 Skip: 0mm
 Matrix: 512 x 192
 FOV: 14cm
 Bandwidth: +/- 16 kHz
 Readout direction: S/I
 Options: Fat Saturation, EDR, ZIP512

Name: Coronal FS 3D SPGR

Pulse Sequence: 3D SPGR
 Scan plane: coronal
 Slices: 32
 TR: 60msec
 TE: 5msec
 Flip angle: 40°
 TI: N/A
 Echo train length: N/A
 NEX: 1
 Slice thickness: 3mm
 Skip: 0mm
 Matrix: 512 x 192
 FOV: 14cm
 Bandwidth: +/- 16 kHz
 Readout direction: S/I
 Options: Fat Saturation, EDR, ZIP512

Name: Sagittal 2D IW (ACL)

Pulse Sequence: 2D FSE
 Scan plane: sagittal oblique (parallel to ACL)
 Slices: 18 (coverage should include bony connections)
 TR: 2000msec
 TE: 40msec
 Flip angle: N/A
 TI: N/A
 Echo train length: 8
 NEX: 2
 Slice thickness: 2mm, interleaved
 Skip: 0mm
 Matrix: 256 x 192
 FOV: 16cm
 Bandwidth: +/- 20-21 kHz
 Readout direction: S/I
 Options: NPW, VBw, EDR, TRF, Fast, ZIP512

Siemens 1.5T MR Exam Protocol

Name: 3-plane Localizer

Pulse Sequence:	2D FL
Scan plane:	axial, coronal, sagittal
Slices:	5 coronal, 5 sagittal, 7 axial (17 total)
TR:	60msec
TE:	10msec
Flip angle:	50°
TI:	N/A
Echo train length:	N/A
NEX:	1
Slice thickness:	5mm
Distance Factor:	20%
Matrix:	256 x 128
FOV:	24cm
Bandwidth:	150 Hz/pixel
Phase direction:	Sagittal A/P, Coronal R/L, Axial A/P
Options:	none

Name: Sagittal 2D IW (ACL)

Pulse Sequence:	2D TSE
Scan	plane:sagittal oblique (parallel to ACL)
Slices:	18 (coverage should include bony connections)
TR:	2000msec
TE:	40msec
Flip angle:	N/A
TI:	N/A
Echo train length:	8
NEX:	2
Slice thickness:	2mm, interleaved
Distance Factor:	0%
Matrix:	256 x 192
FOV:	16cm
Bandwidth:	100% Hz/pixel
Phase direction:	A/P
Options:	100% phase oversampling, elliptical filter, Large FOV filter

Name: Coronal 2D PD FS

Routine

Slice group	1
Slices	28
Dist. factor	15 [%]
Position	-
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0 [deg]
Phase oversampling	100 [%]
FoV read	140 [mm]
FoV phase	100.0 [%]
Slice thickness	3 [mm]
TR	4500 [ms]
TE	Minimum full TE (~15ms)
Averages	1
Concatenations	1
Filter	Elliptical Filter
Coil elements	EX

Contrast

MTC	Off
Magn. preparation	None
Flip angle	180 [deg]
Reconstruction	Magnitude
Fat suppr.	Fat sat
Fat sat. mode	Strong
Water suppr.	None
Measurements	1

Resolution

Base resolution	512
Phase resolution	50 [%]
Phase partial Fourier	Off
Filter 1	
Raw filter	Off
Filter 2	
Large FoV	On
Filter 3	
Normalize	Off
Filter 4	
Elliptical filter	On
Interpolation	Off

PAT mode	None

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.	None

System

Save uncombined	Off
Scan at current TP	Off
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	P >> A
Transversal	H >> F
CP Extremity / EX	On
Body	Off

Shim mode	Standard
Confirm freq. adjustment	Off
Assume Silicone	Off
Ref. amplitude [1H]	180.000 [V]
Adjust volume	
Position	
Orientation	Coronal
Rotation	0 [deg]
F >> H	140 [mm]
R >> L	140 [mm]
A >> P	97 [mm]

Physio

1st Signal/Mode	None

Dark blood	Off

Resp. control	Off

Inline

Subtract	0
Std-Dev-Sag	0
Std-Dev-Cor	0
Std-Dev-Tra	0
Std-Dev-Time	0
MIP-Sag	0
MIP-Cor	0
MIP-Tra	0
MIP-Time	0
Save original images	1

Sequence

Introduction	On
Dimension	2D
Compensate T2 decay	Off
Averaging mode	Long term
Contrasts	1
Bandwidth	100 [Hz/Px]
Flow comp.	No
Allowed delay	30 [s]
Echo spacing	16.4 [ms]

Turbo factor	8
RF pulse type	Fast
Gradient mode	Normal

Name: Sagittal 2D PD FS

Routine

Slice group	1
Slices	28
Dist. factor	15 [%]
Position	-
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0 [deg]
Phase oversampling	100 [%]
FoV read	140 [mm]
FoV phase	100.0 [%]
Slice thickness	3 [mm]
TR	4500 [ms]
TE	Minimum full TE (~15ms)
Averages	1
Concatenations	1
Filter	Elliptical Filter
Coil elements	EX

Contrast

MTC	Off
Magn. preparation	None
Flip angle	180 [deg]
Reconstruction	Magnitude
Fat suppr.	Fat sat
Fat sat. mode	Strong
Water suppr.	None
Measurements	1

Resolution

Base resolution	512
Phase resolution	50 [%]
Phase partial Fourier	Off
Filter 1	
Raw filter	Off
Filter 2	
Large FoV	On
Filter 3	
Normalize	Off
Filter 4	

Elliptical filter	On
Interpolation	Off

PAT mode	None

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.	None

System

Save uncombined	Off
Scan at current TP	Off
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	P >> A
Transversal	H >> F
CP Extremity / EX	On
Body	Off

Shim mode	Standard
Confirm freq. adjustment	Off
Assume Silicone	Off
Ref. amplitude [1H]	180.000 [V]
Adjust volume	
Position	-
Orientation	Sagittal
Rotation	0 [deg]
F >> H	140 [mm]
R >> L	140 [mm]
A >> P	97 [mm]

Physio

1 st Signal/Mode	None

Dark blood	Off

Resp. control	Off

Inline

Subtract	0
Std-Dev-Sag	0
Std-Dev-Cor	0
Std-Dev-Tra	0
Std-Dev-Time	0
MIP-Sag	0
MIP-Cor	0
MIP-Tra	0
MIP-Time	0
Save original images	1

Sequence

Introduction	On
Dimension	2D
Compensate T2 decay	Off
Averaging mode	Long term
Contrasts	1
Bandwidth	100 [Hz/Px]
Flow comp.	No
Allowed delay	30 [s]
Echo spacing	16.4 [ms]

Turbo factor	8
RF pulse type	Fast
Gradient mode	Normal

Name: Sagittal 3D SPGR FS

Routine

Slab group	1
Slabs	1
Dist. factor	20 [%]
Position	-
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0 [deg]
Phase oversampling	20 [%]
Slice oversampling	20 [%]
Slices per slab	32
FoV read	140 [mm]
FoV phase	100 [%]
Slice thickness	3 [mm]
TR	60 [ms]
TE	5.58 [ms]
Averages	1
Concatenations	1
Filter	Elliptical Filter
Coil elements	EX

Contrast

MTC	Off
Magn. preparation	None
Flip angle	40 [deg]
Reconstruction	Magnitude
Fat suppr.	Fat sat (or Water Excitation)
Water suppr.	None
Measurements	1

Resolution

Base resolution	512
Phase resolution	37.5 [%]
Slice resolution	100 [%]
Phase partial Fourier	8/8
Slice partial Fourier	8/8
Filter 1	
Raw filter	Off
Filter 2	
Large FoV	On
Filter 3	
Normalize	Off
Filter 4	
Elliptical filter	On
Interpolation	Off

PAT mode	None

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Saturation mode	Standard
Special sat.	None

System

Save uncombined	Off
Scan at current TP	Off
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	P >> A
Transversal	H >> F
CP Extremity / EX	On
Body	Off

Shim mode	Standard
Confirm freq. adjustment	Off
Assume Silicone	Off
Ref. amplitude [1H]	180.000 [V]
Adjust volume	
Position	-
Orientation	Sagittal
Rotation	0 [deg]
F >> H	140 [mm]
R >> L	140 [mm]
A >> P	96 [mm]

Physio

1 st Signal/Mode	None
Segments	1

Guidelines	None
Dark blood	Off

Resp. control	Off

Inline

Subtract	0
Std-Dev-Sag	0
Std-Dev-Cor	0
Std-Dev-Tra	0
Std-Dev-Time	0
MIP-Sag	0
MIP-Cor	0
MIP-Tra	0
MIP-Time	0
Save original images	1

Wash - In	0
Wash - Out	0
TTP	0
PEI	0
MIP - time	0

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	On
Phase stabilisation	Off
Averaging mode	Short term
Asymmetric echo	Not allowed
Contrasts	1
Bandwidth	130 [Hz/Px]
Flow comp.	No

RF pulse type	Fast
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On

Name: Coronal 3D SPGR FS

Routine

Slab group	1
Slabs	1
Dist. factor	20 [%]
Position	-
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0 [deg]
Phase oversampling	20 [%]
Slice oversampling	20 [%]
Slices per slab	32
FoV read	140 [mm]
FoV phase	100 [%]
Slice thickness	3 [mm]
TR	60 [ms]
TE	5.58 [ms]
Averages	1
Concatenations	1
Filter	Elliptical Filter
Coil elements	EX

Contrast

MTC	Off
Magn. preparation	None
Flip angle	40 [deg]
Reconstruction	Magnitude
Fat suppr.	Fat sat (or Water Excitation)
Water suppr.	None
Measurements	1

Resolution

Base resolution	512
Phase resolution	37.5 [%]
Slice resolution	100 [%]
Phase partial Fourier	8/8
Slice partial Fourier	8/8
Filter 1	
Raw filter	Off
Filter 2	
Large FoV	On
Filter 3	
Normalize	Off
Filter 4	
Elliptical filter	On
Interpolation	Off

PAT mode	None

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Saturation mode	Standard
Special sat.	None

System

Save uncombined	Off
Scan at current TP	Off
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T

Sagittal	R >> L
Coronal	P >> A
Transversal	H >> F
CP Extremity / EX	On
Body	Off

Shim mode	Standard
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	180.000 [V]
Adjust volume	
Position	-
Orientation	Coronal
Rotation	0 [deg]
F >> H	140 [mm]
R >> L	140 [mm]
A >> P	96 [mm]

Physio

1 st Signal/Mode	None
Segments	1

Guidelines	None
Dark blood	Off

Resp. control	Off

Inline

Subtract	0
Std-Dev-Sag	0
Std-Dev-Cor	0
Std-Dev-Tra	0
Std-Dev-Time	0
MIP-Sag	0
MIP-Cor	0
MIP-Tra	0
MIP-Time	0
Save original images	1

Wash - In	0
Wash - Out	0
TTP	0
PEI	0
MIP - time	0

Sequence

Introduction	1
Dimension	3D
Elliptical scanning	On
Phase stabilisation	Off
Averaging mode	Short term
Asymmetric echo	Not allowed
Contrasts	1
Bandwidth	130 [Hz/Px]
Flow comp.	No

RF pulse type	Fast
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On

Philips 1.5T MR Exam Protocol

Name: 3-plane Localizer at RF coil isocenter (R/L60, A0)

Pulse Sequence: ffe/m, 2D, gr
 Scan plane: axial, coronal, sagittal
 Slices: 5 coronal, 5 sagittal,
 7 axial (17 total)
 TR: 9msec
 TE: 4.6msec
 Flip angle: 50°
 TI: N/A
 Echo train length: N/A
 NEX: 1
 Slice thickness: 5mm
 Skip: 6mm
 Matrix: 256 x 128
 (50% phase sampling)
 FOV: 24cm
 Bandwidth: 150 Hz/pixel
 Phase direction: Sagittal A/P,
 Coronal R/L, Axial A/P
 Options: none

TSE factor = 8
 startup echoes = 0
 profile order = linear
 DRIVE = no
 ultrashort = yes
 shift = 0.00
 Echoes = 1
 partial echo = no
 TE = min full
 Flip angle (deg) = 90.00
 Refocusing control = yes
 TR = 4500.0
 Halfscan = no
 Water-fat shift = user defined
 (pixels) = 2.20
 Shim = standard
 SPIR = no
 SPAIR = no
 BB pulse = yes
 (??? *Bionomial Fat Sat*)

Name: Sagittal 2D PD FS

Nucleus = H1
 Coil selection = SENSE-Knee-8
 ch. combination = 12345678
 connection = d
 Dual coil = no
 Homogeneity correction = none
 CLEAR = no
 Nucleus = H1
 FOV FH (mm) = 140
 AP (mm) = 140
 Matrix = 512
 Scan percentage = 50.0
 Slice thickness (mm) = 3.00
 Fold-over suppression = yes
 % = 70
 Reconstruction matrix = 512
 SENSE = no
 Stacks = 1
 type = parallel
 Slices = 28
 slice gap = user defined
 gap (mm) = 0.0
 Slice orientation = sagittal
 Fold-over direction = AP
 Fat shift direction = F
 Minimum n° of packages = 1
 Slice scan order = default
 PlanAlign = no
 REST slabs = 0
 Catheter tracking = no
 Interactive positioning = no
 Allow table movement = no
 Patient position = feet first
 orientation = supine
 Scan type = Imaging
 Scan mode = MS
 technique = SE
 Modified SE = no
 Acquisition mode = cartesian
 Fast Imaging mode = TSE
 shot mode = multishot

ProSet = no
 MTC = no
 Research prepulse = no
 Zoom imaging = no
 Diffusion mode = no
 Elastography mode = no
 SAR mode = high
 B1 mode = user defined
 amplitude (uT) = 13.00
 (may need to lower to 11.0)
 PNS mode = moderate
 Gradient mode = maximum
 SofTone mode = no
 Cardiac synchronization = no
 Respiratory compensation = no
 Navigator respiratory comp = no
 Flow compensation = no
 Temporal slice spacing = default
 Motion smoothing = no
 NSA = 1
 Manual start = no
 Dynamic study = no
 Flow labelling = none
 Preparation phases = auto
 Gradient demo = no
 SENSE ref. scan = no
 SmartPlan scout = no
 B0 field map = no
 B1 field map = no
 MIP/MPR = no
 Images = M
 Autoview image = M
 Calculated images = no
 Reference tissue = Skeletal muscle
 Preset window contrast = soft
 Reconstruction mode = immediate
 Save raw data = no
 Push to workstation = no
 Hardcopy protocol = no
 Ringing filtering = no
 Geometry correction = default

Name: Coronal 2D PD FS

Analogous to Sagittal 2D PD FS, except with:

Slice orientation = coronal
 Fold-over direction = RL

Shim = standard
 SPIR = no
 SPAIR = no
 BB pulse = yes
 (??? Binomial Fat Sat)

Name: Sagittal 3D SPGR FS

Nucleus = H1
 Coil selection = SENSE-Knee-8
 ch. combination = 12345678
 connection = d
 Dual coil = no
 Homogeneity correction = none
 CLEAR = no
 Nucleus = H1
 FOV FH (mm) = 140
 AP (mm) = 140
 Matrix = 512
 Scan percentage = 37.5
 Slice thickness (mm) = 3.00
 Fold-over suppression = no
 Reconstruction matrix = 512
 SENSE = no
 Stacks = 1
 type = parallel
 Slices = 32
 slice gap = user defined
 gap (mm) = 0.0
 Slice orientation = sagittal
 Fold-over direction = AP
 Fat shift direction = F
 Minimum n° of packages = 1
 Slice scan order = default
 PlanAlign = no
 REST slabs = 0
 Catheter tracking = no
 Interactive positioning = no
 Allow table movement = no
 Patient position = feet first
 orientation = supine
 Scan type = Imaging
 Scan mode = 3D
 technique = ffe/m
 Modified SE = no
 Acquisition mode = cartesian
 Fast Imaging mode = GR
 shot mode = multishot
 TSE factor = 1
 Echoes = 1
 partial echo = no
 TE = user defined
 multi (ms) = 4.60
 Flip angle (deg) = 40.00
 Refocusing control = yes
 TR = 60.0
 Halfscan = no
 Water-fat shift = user defined
 (pixels) = 1.750

ProSet = no
 MTC = no
 Research prepulse = no
 Zoom imaging = no
 Diffusion mode = no
 Elastography mode = no
 SAR mode = high
 B1 mode = user defined
 amplitude (uT) = 13.00
 PNS mode = moderate
 Gradient mode = maximum
 SofTone mode = no
 Cardiac synchronization = no
 Respiratory compensation = no
 Navigator respiratory comp = no
 Flow compensation = no
 Temporal slice spacing = default
 Motion smoothing = no
 NSA = 1
 Manual start = no
 Dynamic study = no
 Flow labelling = none
 Preparation phases = auto
 Gradient demo = no
 SENSE ref. scan = no
 SmartPlan scout = no
 B0 field map = no
 B1 field map = no
 MIP/MPR = no
 Images = M
 Autoview image = M
 Calculated images = no
 Reference tissue = Skeletal muscle
 Preset window contrast = soft
 Reconstruction mode = immediate
 Save raw data = no
 Push to workstation = no
 Hardcopy protocol = no
 Ringing filtering = no
 Geometry correction = default

Name: Coronal 3D SPGR FS

Analogous to Sagittal 3D SPGR FS, except with

Slice orientation = coronal
 Fold-over direction = RL

Name: Sagittal 2D IW (ACL)

Nucleus	= H1	Water-fat shift	= user defined
Coil selection	= SENSE-Knee-8	(pixels)	= 2.20
ch. combination	= 12345678	(if fuzzy, decrease to 1.8)	
connection	= d	Shim	= standard
Dual coil	= no	SPIR	= no
Homogeneity correction	= none	SPAIR	= no
CLEAR	= no	BB pulse	= no
Nucleus	= H1	ProSet	= no
FOV		MTC	= no
FH (mm)	= 140	Research prepulse	= no
AP (mm)	= 140	Zoom imaging	= no
Matrix	= 256	Diffusion mode	= no
Scan percentage	= 75.0	Elastography mode	= no
Slice thickness (mm)	= 2.00	SAR mode	= high
Fold-over suppression	= yes	B1 mode	= user defined
%	= 70	amplitude (uT)	= 13.00
Reconstruction matrix	= 512	(may need to lower to 11.0)	
SENSE	= no	PNS mode	= moderate
Stacks	= 1	Gradient mode	= maximum
type	= parallel	Softone mode	= no
Slices	= 18	Cardiac synchronization	= no
(enough to cover ACL insertions)		Respiratory compensation	= no
slice gap	= user defined	Navigator respiratory comp	= no
gap (mm)	= 0.0	Flow compensation	= no
Slice orientation	= sagittal oblique	Temporal slice spacing	= default
(parallel to ACL)		Motion smoothing	= no
Fold-over direction	= AP	NSA	= 1
Fat shift direction	= F	Manual start	= no
Minimum n° of packages	= 1	Dynamic study	= no
Slice scan order	= default	Flow labelling	= none
PlanAlign	= no	Preparation phases	= auto
REST slabs	= 0	Gradient demo	= no
Catheter tracking	= no	SENSE ref. scan	= no
Interactive positioning	= no	SmartPlan scout	= no
Allow table movement	= no	B0 field map	= no
Patient position	= feet first	B1 field map	= no
orientation	= supine	MIP/MPR	= no
Scan type	= Imaging	Images	= M
Scan mode	= MS	Autoview image	= M
technique	= SE	Calculated images	= no
Modified SE	= no	Reference tissue	= Skeletal muscle
Acquisition mode	= cartesian	Preset window contrast	= soft
Fast Imaging mode	= TSE	Reconstruction mode	= immediate
shot mode	= multishot	Save raw data	= no
TSE factor	= 8	Push to workstation	= no
startup echoes	= 0	Hardcopy protocol	= no
profile order	= linear	Ringing filtering	= no
DRIVE	= no	Geometry correction	= default
ultrashort	= yes		
shift	= 0.00		
Echoes	= 1		
partial echo	= no		
TE	= user defined		
multi (ms)	= 40.00		
Flip angle (deg)	= 90.00		
Refocusing control	= yes		
TR	= 2000.0		
Halfscan	= no		

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